

nanoLiter Cool Wave™ Dispenser for MS Sample Intro/Prep; Fluidics R&D & More! ONE DEVICE FOR BOTH TECHNIQUES! ESI MALDI

NEW 21st century sample handling, prep fluidic tool set for ESI, MALDI, Sample Prep and R&D and other simple fluidic apps.

nLs yields *major* sensitivity, reproducibility increases for ESI,MALDI with the SAME device! Dispense *viscous* liquids, cells, blood, serum, polymers into/near any ESI! nanoLiter Cool Wave is the lowest cost, fastest, most efficient and versatile, ESI sample introduction, MALDI dispensing tool in the world, literally! Shoot droplets in a line like this..........! Get it?




G. Groenewold, et al, INL, Analytical, Chem, 6/13.



IBF shoots sample into any ESI MS. Also append IBF to ANY robotic dispenser.






1. 0 msec!



Ca. 35 nLs




J. Harmon et al, USF

Fastest, most efficient MS sample introduction/prep. in world, literally! Makes excellent MALDI, SIMS, LDI crystals. See M. Gross, et al JASMS 08/08. Morph your **pump** into an ion source !!! that can also be a non-touch dispenser!!! Shoot cells, liquids into your MS. Do droplet research. Is ESI 100x more sensitive? 100% IBF ESI/LC/MS is more sensitivity than nano-Spray! ASMS 2014, in press.

Save precious samples, expensive labeled standards. Make experiments safer, produce 10-100x less waste. IBF is GREEN tech! With full purchase, we often consult and train.

License IP, including international positions available.



MS syringes with 100 removable tips, \$149.00. Ask about Pipette/MS. Explore LO 3D printing.